



Recombinant *Saccharomyces cerevisiae* Stationary phase gene 1 protein (SPG1)

Product Code	CSB-CF413858STA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A6ZUP4
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Sequence	MKLDSGIYSEAQRVVRTPKFRYIMLGLVGAAVVPTAYMRRGYTVPAHSLDNIN GVDTTKA SVMGTEQRAAMTKGKSLQEMMDDDEVTYLMFSSIM
Source	in vitro <i>E.coli</i> expression system
Target Names	SPG1
Protein Names	Recommended name: Stationary phase gene 1 protein
Expression Region	1-95
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.