



Recombinant Na (+)/H (+) antiporter subunit F1

Product Code	CSB-CF352136FKZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P60694
Product Type	Transmembrane Protein
Immunogen Species	Staphylococcus aureus
Sequence	MNHNVIIVIALIIVVISMLAMLIRVVLGPSLADRVALDAIGLQLMAVIALFSILLNIK Y MIVVIMMIGILAF LGTAVFSKFM DKGKVIEHDQNHTD
Source	in vitro E.coli expression system
Target Names	mnhF1
Protein Names	Recommended name: Na(+)/H(+) antiporter subunit F1 Alternative name(s): Mnh complex subunit F1 Mrp complex subunit F1
Expression Region	1-97
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.