



Recombinant Gorilla gorilla gorilla NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 3 (NDUFB3)

Product Code	CSB-CF610337GGZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0MQD1
Product Type	Transmembrane Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Sequence	AHEHGHEHGHHKMELPDYRQWKIEGTPLETIQKKLAAKGLRDPWGRNEAWR YMGGFAKSV SFSDVFFKGFKWGFAAFVAVGAEYYLESLNKDKKHH
Source	in vitro E.coli expression system
Target Names	NDUFB3
Protein Names	Recommended name: NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 3 Alternative name(s): Complex I-B12 Short name= CI-B12 NADH-ubiquinone oxidoreductase B12 subunit
Expression Region	2-98
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.