



Recombinant *Platyrrhinus helleri* NADH-ubiquinone oxidoreductase chain 4L (MT-ND4L)

Product Code	CSB-CF624918PBAI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q1HV02
Product Type	Transmembrane Protein
Immunogen Species	<i>Platyrrhinus helleri</i> (Heller's broad-nosed bat) (<i>Vampyrops helleri</i>)
Sequence	MSLTYMNMFMFAFTISLLGLLLYRSHMMSSLLCLEGMMLSLFVMMTMVILNTHL TLASMIP IILLVFAACEAALGLSLLVMVSTTYGMDYVQNLNLLQC
Source	in vitro E.coli expression system
Target Names	MT-ND4L
Protein Names	Recommended name: NADH-ubiquinone oxidoreductase chain 4L EC= 1.6.5.3 Alternative name(s): NADH dehydrogenase subunit 4L
Expression Region	1-98
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.