



Recombinant *Sminthopsis crassicaudata* NADH-ubiquinone oxidoreductase chain 4L (MT-ND4L)

Product Code	CSB-CF657013SZT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q32US1
Product Type	Transmembrane Protein
Immunogen Species	<i>Sminthopsis crassicaudata</i> (Fat-tailed dunnart) (<i>Phascogale crassicaudata</i>)
Sequence	MLSINLNLIVAFLLALMGVLIYRSHLMSTLLCLEGMMLSLFILMTLLITHFHMFSM SMTP LILLVFSACEAAIGLALLVKISATHGSDHIQNLNLLQC
Source	in vitro <i>E.coli</i> expression system
Target Names	MT-ND4L
Protein Names	Recommended name: NADH-ubiquinone oxidoreductase chain 4L EC= 1.6.5.3 Alternative name(s): NADH dehydrogenase subunit 4L
Expression Region	1-98
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.