




Recombinant Human Ribonuclease kappa (RNASEK)

Product Code	CSB-CF761293HU
Abbreviation	Recombinant Human RNASEK protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6P5S7
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose
Product Type	Transmembrane Proteins
Immunogen Species	Homo sapiens (Human)
Sequence	MASLLCCGPKLAACGIVLSAWGVIMLIMLGIFFNVHSAVLIEDVPFTEKDFENG PQNIYN LYEQVSYNCFIAAGLYLLLGGFSFCQVRLNKRKEYMVR
Source	in vitro E.coli expression system
Target Names	RNASEK
Protein Names	Recommended name: Ribonuclease kappa Short name= RNase K Short name= RNase kappa EC= 3.1.-.-
Expression Region	1-98
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Protein Length	full length protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.