



Recombinant *Macaca fascicularis* NADH-ubiquinone oxidoreductase chain 4L (MT-ND4L)

Product Code	CSB-CF765947MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q7J3C3
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sequence	MTPTYMNIMLAFTISLLGMLTYRSHLVASLLCLEGMMMSLFIMATLIASNTHFPL VNIMP IILLVFAACETAVGLALLISISNTYGLDYVHNLNLLQC
Source	in vitro <i>E.coli</i> expression system
Target Names	MT-ND4L
Protein Names	Recommended name: NADH-ubiquinone oxidoreductase chain 4L EC= 1.6.5.3 Alternative name(s): NADH dehydrogenase subunit 4L
Expression Region	1-98
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.