



Recombinant Drosophila melanogaster Signal peptidase complex subunit 1 (Spase12)

Product Code	CSB-CF893724DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VAL0
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MLDIQTHMDFAGQGKAERWSRFITFFGIVGLVYGAFVQQFSQTVYILGAGFVL SSLITI PPWPLYRRNALKWQKPIDTDAKSSSESSESGDEGKKKKKQ
Source	in vitro E.coli expression system
Target Names	Spase12
Protein Names	Recommended name: Signal peptidase complex subunit 1 EC= 3.4.-.- Alternative name(s): Microsomal signal peptidase 12 kDa subunit Short name= SPase 12 kDa subunit
Expression Region	1-98
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.