



Recombinant *Deinococcus deserti* NADH-quinone oxidoreductase subunit K (nuoK)

Product Code	CSB-CF509018DIY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C1D0I1
Product Type	Transmembrane Protein
Immunogen Species	<i>Deinococcus deserti</i> (strain VCD115 / DSM 17065 / LMG 22923)
Sequence	MVPTTYYLALSGLLFALGMIGVLTRRTAIMVFLSVELMLNAANLSLVAFARAWG DLTGQT AVFIVMTLAAA EVAIGLAIIVAIFRKRETTNVDLGLKG
Source	in vitro E.coli expression system
Target Names	nuoK
Protein Names	Recommended name: NADH-quinone oxidoreductase subunit K EC= 1.6.99.5 Alternative name(s): NADH dehydrogenase I subunit K NDH-1 subunit K
Expression Region	1-100
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.