



Recombinant Chicken Protein Asterix (WDR83OS)

Product Code	CSB-CF517258CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	F8RT80
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	MADPRRPARVTRYKPPTTESNPALEDPTPDYMNLLGMVFSMCGMLMLKWKW AWIAVYCSF ISFANSRSEDTKQMMSSFMLSISAVVMSYLNPNPQPMSPPW
Source	in vitro E.coli expression system
Target Names	WDR83OS
Protein Names	Recommended name: Protein Asterix Alternative name(s): Protein WDR83OS homolog
Expression Region	1-101
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.