



Recombinant Rat Gamma-secretase subunit PEN-2 (Psenen)

Product Code	CSB-CF757686RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6QI68
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MNLERSVSNEEKLNLCKRKYLLGGFAFLPFLWLVNIFWFFKEAFFAPAYTEQSQIK GYVWRS AVGFLLFWIVLTTWITIFQIYRPRWGALGDYLSFTIPLGTP
Source	in vitro E.coli expression system
Target Names	Psenen
Protein Names	Recommended name: Gamma-secretase subunit PEN-2 Alternative name(s): Liver regeneration-related protein LRRGT00140 Presenilin enhancer protein 2
Expression Region	1-101
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.