



Recombinant Maize streak virus genotype E Movement protein (V2)

Product Code	CSB-CF835554MIAM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q91MG4
Product Type	Transmembrane Protein
Immunogen Species	Maize streak virus genotype E (isolate Pat) (MSV)
Sequence	MDPQSAVYSLPRVPTAAPPNAGVPWSHVGEVAVLSFVALICIIYLLYLWVLRDLI LVLKAR RGRSTEELIFGSEAVDRRSPIPTLEPTAPVHPGPFVPGSG
Source	in vitro E.coli expression system
Target Names	V2
Protein Names	Recommended name: Movement protein Short name= MP
Expression Region	1-101
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.