



Recombinant Mouse Gamma-secretase subunit PEN-2 (Psenen)

Product Code	CSB-CF861450MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9CQR7
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MNLERSVSNEEKLNLCKRKYLLGGFAFLPFLWLVNIFWFFREAFAPAYTEQSQIK GYVWRS AVGFLLFWIILATWITIFQIYRPRW GALGDYLSFTIPLGTP
Source	in vitro E.coli expression system
Target Names	Psenen
Protein Names	Recommended name: Gamma-secretase subunit PEN-2 Alternative name(s): Presenilin enhancer protein 2
Expression Region	1-101
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.