



Recombinant Schizosaccharomyces pombe UPF0136 membrane protein P14E8.05c (SPAP14E8.05c)

Product Code	CSB-CF868443SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9P7G3
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MTPDQNALVLSFLLTVGGLIGYLRKKSIVSLIAGTALGANFAWASKLMERGSS QGINYAF YGSLVLLASSGPRFYKSRKPVPMILTVLGVISTWYFYRLWA
Source	in vitro E.coli expression system
Target Names	SPAP14E8.05c
Protein Names	Recommended name: UPF0136 membrane protein P14E8.05c
Expression Region	1-101
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.