



# Recombinant Nostoc sp. NAD (P)H-quinone oxidoreductase subunit 4L (ndhE)

<b>Product Code</b>	CSB-CF894848NHR
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q9WWM4
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576)
<b>Sequence</b>	MQLRYFLLLAALFCIGIYGLITSRNAVRLMSIELLLNAVNLNLMAFSNYLDSTL IKGQ VFTVFVITVAAAEEAAVGLAIVLAIYRNRDVTVDMEQFNLLKW
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	ndhE
<b>Protein Names</b>	Recommended name: NAD(P)H-quinone oxidoreductase subunit 4L EC= 1.6.5.- Alternative name(s): NAD(P)H dehydrogenase subunit 4L NADH-plastoquinone oxidoreductase subunit 4L NDH-1, subunit 4L NDH-E
<b>Expression Region</b>	1-101
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.