



Recombinant Human Protein PAPPAS (PAPPA-AS1)

Product Code	CSB-CF719029HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5QFB9
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MRYGFVRKKHRGLFLTVAALPIWNPISSEFVKWYKSHKLSQHCIRICGHLCCQK HLDMFLS VIGQRWPIDVFSSVFDHQVSAIGSDIWWFLKFLVSSFFFF
Source	in vitro E.coli expression system
Target Names	PAPPA-AS1
Protein Names	Recommended name: Protein PAPPAS Alternative name(s): DIPLA1 antisense RNA 1 DIPLA1 antisense gene protein 1 DIPLA1-antisense expressed PAPPAS antisense RNA 1 PAPPAS antisense gene protein 1 PAPPAS-antisense expressed
Expression Region	1-102
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.