



Recombinant Pig Signal peptidase complex subunit 1 (SPCS1)

Product Code	CSB-CF022514PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B0FWK4
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MLEHLSSLPTQMDYKGQKLAEQMFQGII LFS AIVGFIYGYVAEQFGWTVYIVMA GFAFSC LLTLPPWPIYRRHPLKWLPVQDSGTEDKKPGERKIKRHAKNN
Source	in vitro E.coli expression system
Target Names	SPCS1
Protein Names	Recommended name: Signal peptidase complex subunit 1 EC= 3.4.-.- Alternative name(s): Microsomal signal peptidase 12 kDa subunit Short name= SPase 12 kDa subunit
Expression Region	1-102
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.