



Recombinant Human Succinate dehydrogenase [ubiquinone] cytochrome b small subunit, mitochondrial (SDHD)

Product Code	CSB-CF020909HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14521
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	SGSKAASLHWTSERVVSVLLLGLLPAAYLNPCSAMDYSLAAALTLHGHWGLG QVVTDYVH GDALQKAAKAGLLALSALTFAGLCYFNYHDVGICKAVAMLWKL
Source	in vitro E.coli expression system
Target Names	SDHD
Protein Names	Recommended name: Succinate dehydrogenase [ubiquinone] cytochrome b small subunit, mitochondrial Short name= CybS Alternative name(s): CII-4 QPs3 Succinate dehydrogenase complex subunit D Succinate-ubiquinone oxidoreductase cyto
Expression Region	57-159
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.