



Recombinant Pan troglodytes Glycophorin-B (GYPB)

Product Code	CSB-CF631735EQV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q28914
Product Type	Transmembrane Protein
Immunogen Species	Pan troglodytes (Chimpanzee)
Sequence	SSTTEVAMHTSTSSSVTKSYISSQTNDKHKGDTYPATLGAHEVSEISVTTVYPP EEDNGE WVQPVHPFSRPAPVVIIILCVMAGVIGTILLISYGIRLLIKA
Source	in vitro E.coli expression system
Target Names	GYPB
Protein Names	Recommended name: Glycophorin-B Alternative name(s): CD_antigen= CD235b
Expression Region	20-123
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.