



Recombinant Danio rerio UPF0697 protein C8orf40 homolog (zgc:114040)

Product Code	CSB-CF687180DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q4V921
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MGGYGVMADDEESLDYSVHEAWNEATNVYLLVILVSFALLMYARKNKRKIMRIF TLPPTVG SSSEPNFYDSLQKVRLRQQLEMYSLARKYDQQQSQSESQVQLSME
Source	in vitro E.coli expression system
Target Names	smim19
Protein Names	Recommended name: UPF0697 protein C8orf40 homolog
Expression Region	1-104
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.