



Recombinant Mitochondrial inner membrane organizing system protein F54A3.5 (F54A3.5)

Product Code	CSB-CF885645CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9N4K0
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MAPGASAPAASRSEDEVGQKIDRCFADSLKVTGGVAIGIVASVAFFKSRSWPI WFGSGV GLGTGWSNCRHDFASPYVLHGKRVPAQDSQGKPAYNIITEQHKQ
Source	in vitro E.coli expression system
Target Names	F54A3.5
Protein Names	Recommended name: Mitochondrial inner membrane organizing system protein F54A3.5
Expression Region	1-105
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.