



Recombinant *Drosophila yakuba* Vacuolar ATPase assembly integral membrane protein VMA21 (GE19686)

Product Code	CSB-CF025866DMR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B4PET6
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila yakuba</i> (Fruit fly)
Sequence	MSTKNKKAAGGNGGAPKQTRQQSHDSQDYSSFKTVLFYCMLIVFLPVLTFVVLKGFVLDQ FLDISEVKVNIASAVGAVVALHVALGLYIYRAYFGAPGSKASKTD
Source	in vitro <i>E.coli</i> expression system
Target Names	GE19686
Protein Names	Recommended name: Vacuolar ATPase assembly integral membrane protein VMA21
Expression Region	1-105
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.