



Recombinant *Saccharomyces cerevisiae* Uncharacterized membrane protein YIL054W (YIL054W)

Product Code	CSB-CF340267SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40524
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MAPKAFFVCLPWVLPVPRHALIVRQAGNPYHFLAYTNPRAPGKLQDSHCPVFFM GIIIIITII TVTLAIIIIINIIFLTLFDDGMCFYCSLLTFSFVSFNFDHFDHFDL
Source	in vitro <i>E.coli</i> expression system
Target Names	YIL054W
Protein Names	Recommended name: Uncharacterized membrane protein YIL054W
Expression Region	1-105
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.