



Recombinant *Saccharomyces cerevisiae* Altered inheritance of mitochondria protein 5, mitochondrial (AIM5)

Product Code	CSB-CF511318SVL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C7GMW8
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain JAY291) (Baker's yeast)
Sequence	MSKLGPLARSVKWTLSVGVIGSVFYLYRYSNNGYFYDHDATWLKQDHQVQDL VDRKEVVP GETRNRKLVVTTDDGTAWSRTMGESIKDIWNEQIRNSVDWIYSWGKN
Source	in vitro E.coli expression system
Target Names	AIM5
Protein Names	Recommended name: Altered inheritance of mitochondria protein 5, mitochondrial Alternative name(s): Found in mitochondrial proteome protein 51
Expression Region	1-106
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.