



Recombinant Mouse PAT complex subunit Asterix (Wdr83os)

Product Code	CSB-CF744431MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZWX0
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Transmembrane Proteins
Immunogen Species	Mus musculus (Mouse)
Sequence	MSTNNMSDPRRPNKVLRYPKPPPSECNPALDDPTPDYMNLLGMIFSMCGLMLK LKWCAWVA VYCSFISFANSRSEDTKQMMSSFMLSISAVVMSYLQNPQPMTTPPW
Source	in vitro E.coli expression system
Target Names	Wdr83os
Protein Names	Recommended name: Protein Asterix Alternative name(s): Protein WDR83OS homolog
Expression Region	1-106
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.