



Recombinant *Saccharomyces cerevisiae* Putative uncharacterized membrane protein YAR047C (YAR047C)

Product Code	CSB-CF331040SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P39557
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MYQTSPLSLFYFQVLVLPKFLECFPCFPYHKISLVALLSFFYCYQLQTNMIILLSQIK RFLY RQIMIALKIKAKKFWFIFKYFNVSCDARLNFELFYIFQTYVSVDSK
Source	in vitro <i>E.coli</i> expression system
Target Names	YAR047C
Protein Names	Recommended name: Putative uncharacterized membrane protein YAR047C
Expression Region	1-106
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.