



Recombinant Enterobacteria phage lambda Holin (S)

Product Code	CSB-CF361129ECW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P03705
Product Type	Transmembrane Protein
Immunogen Species	Escherichia phage lambda (Bacteriophage lambda)
Sequence	MKMPEKHDLLAAILAAKEQGIGAILAFAMAYLRGRYNGGAFTKTVIDATMCAIIA WFIRD LLDFAGLSSNLAYITSVFIGYIGTDSIGSLIKRFAAKKAGVEDGRNQ
Source	in vitro E.coli expression system
Target Names	S
Protein Names	Recommended name: Holin Alternative name(s): gpS protein Including the following 2 domains: Lysis protein S Lysis inhibitor
Expression Region	1-107
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.