



Recombinant Bat coronavirus 279/2005 Protein 7a (7a)

Product Code	CSB-CF611307BFB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q0Q470
Product Type	Transmembrane Protein
Immunogen Species	Bat coronavirus 279/2005 (BtCoV) (BtCoV/279/2005)
Sequence	ELYHYQECVIRGTTVLLKEPCPSGTYEGNSPFHPLADNKFALTCISTHFAFACA DGTRHTY QLRARSVSPKLFTRQEEVHQELYSPFLFLIVAALVFIILCFTIKRKTE
Source	in vitro E.coli expression system
Target Names	7a
Protein Names	Recommended name: Protein 7a Alternative name(s): Accessory protein 7a
Expression Region	16-122
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.