



Recombinant Human UPF0697 protein C8orf40 (C8orf40)

Product Code	CSB-CF856922HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96E16
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAGGYGVMGDDGSIDYTVHEAWNEATNVYLIVILVSFGLFMYAKRNKRRIMRI FSVPSTE ETLSEPNFYDTISKIRLRQQLEMYSISRKYDYQQPQNQADSVQLSLE
Source	in vitro E.coli expression system
Target Names	SMIM19
Protein Names	Recommended name: UPF0697 protein C8orf40
Expression Region	1-107
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.