



Recombinant Human UPF0258 protein KIAA1024-like (KIAA1024L)

Product Code	CSB-CF012224HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P59773
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSSVMKNNPLYGDLSLEEAMEERKKNPSWTIEEYDKHSLHTNLSGHLKENPN DLRFWLGD MYTPGFDLLKKEEKQEKHSKFCRMGLILLVVISILVTIVTIITFFT
Source	in vitro E.coli expression system
Target Names	KIAA1024L
Protein Names	Recommended name: UPF0258 protein KIAA1024-like
Expression Region	1-107
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.