



Recombinant Neosartorya fischeri Vacuolar ATPase assembly integral membrane protein VMA21 (vma21)

Product Code	CSB-CF025866NEZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1D7K7
Product Type	Transmembrane Protein
Immunogen Species	Neosartorya fischeri (strain ATCC 1020 / DSM 3700 / CBS 544.65 / FGSC A1164 / JCM 1740 / NRRL 181 / WB 181) (Aspergillus fischerianus)
Sequence	MTSRRSQEKSYAEAAAAPPK EAASSDVTPAVPADVIYKLLGFTAAMVVGPIG MYFITVN SGASSTVAGITAAITANLVLFGYIYVAWLDDREEREAAASKKKEKKAQ
Source	in vitro E.coli expression system
Target Names	vma21
Protein Names	Recommended name: Vacuolar ATPase assembly integral membrane protein VMA21
Expression Region	1-107
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.