



Recombinant Ureaplasma parvum serovar 3 UPF0154 protein UU265 (UU265)

Product Code	CSB-CF889448UAO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9PQM6
Product Type	Transmembrane Protein
Immunogen Species	Ureaplasma parvum serovar 3 (strain ATCC 700970)
Sequence	MNFSSFIFKVSDFKSVIHEASDVVTKADLDNANAHTHSLAVGLGIGIVLFLIAGL IIGY FISMKIMKRQLKKNPPISKDTIRMIYQQVGRKPSESQINEIYNRAVKQK
Source	in vitro E.coli expression system
Target Names	UU265
Protein Names	Recommended name: UPF0154 protein UU265
Expression Region	1-109
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.