



Recombinant Escherichia coli UPF0092 membrane protein YajC (yajC)

Product Code	CSB-CF365062ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0ADZ7
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MSFFISDAVAATGAPAQGSPMSLILMLLVFGLIFYFMILRPQQKRTKEHKKLMD SIAKGD EVLTNGLVGRVTKVAENGYIAIALNDTTEVVIK RDFVAAVLPKGTMKAL
Source	in vitro E.coli expression system
Target Names	yajC
Protein Names	Recommended name: UPF0092 membrane protein YajC
Expression Region	1-110
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.