



Recombinant Escherichia coli Protein-export membrane protein SecG (secG)

Product Code	CSB-CF314852ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0AG99
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MYEALLVVFLIVAIGLVGLIMLQQGKGADMGASFGAGASATLFGSSGSGNFMT RMTALLA TLFFIISLVLGNINSNKTNGKSEWENLSAPAKTEQTQPAAPAKPTSDIPN
Source	in vitro E.coli expression system
Target Names	secG
Protein Names	Recommended name: Protein-export membrane protein SecG Alternative name(s): P12 Preprotein translocase band 1 subunit
Expression Region	1-110
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.