



Recombinant Buchnera aphidicola subsp. Acyrthosiphon pisum UPF0092 membrane protein BU134 (BU134)

Product Code	CSB-CF348123BWZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P57234
Product Type	Transmembrane Protein
Immunogen Species	Buchnera aphidicola subsp. Acyrthosiphon pisum (strain APS) (Acyrthosiphon pisum symbiotic bacterium)
Sequence	MSFFIQNANAVVNGTSESSNSYSLIFMAVIFLLIFYFMLFRPQQKKDKKEHKNLIN SLVQG DEVITTSGLLGRIKKITKNGYILLELNETTEVFQKDFIVSLLPKGTLKSL
Source	in vitro E.coli expression system
Target Names	yajC
Protein Names	Recommended name: UPF0092 membrane protein BU134
Expression Region	1-111
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.