



Recombinant Inner membrane protein ydgC (ydgC)

Product Code	CSB-CF359846EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0ACX1
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MGLVIKAALGALVLLIGVLAKTKNYYIAGLIPLFPTFALIAHYIVASERGIEALRAT II FSMWSIIPYFVYLVSLWYFTGMMRLPAAFVGSVACWGISAWVLIICWIKLH
Source	in vitro E.coli expression system
Target Names	ydgC
Protein Names	Recommended name: Inner membrane protein ydgC
Expression Region	1-111
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.