



Recombinant Escherichia coli O139:H28 Probable 4-amino-4-deoxy-L-arabinose-phosphoundecaprenol flippase subunit ArnE (arnE)

Product Code	CSB-CF419506EJD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A7ZP76
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O139:H28 (strain E24377A / ETEC)
Sequence	MIWLTLVFASLLSVAGQLCQKQATCFVAINKRRKHIVLWLGLALACLGLAMVLW LLVLQN VPVGIAYPMLSLNFWVTAAVKLWHEPVSPRHWCGVAFIIGGIVILGSTV
Source	in vitro E.coli expression system
Target Names	arnE
Protein Names	Recommended name: Probable 4-amino-4-deoxy-L-arabinose-phosphoundecaprenol flippase subunit ArnE Short name= L-Ara4N-phosphoundecaprenol flippase subunit ArnE Alternative name(s): Undecaprenyl phosphate-aminoarabinose flippase subunit ArnE
Expression Region	1-111
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.