



Recombinant Escherichia coli O7:K1 Universal stress protein B (uspB)

Product Code	CSB-CF487214EON
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7NNB6
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O7:K1 (strain IAI39 / ExPEC)
Sequence	MISTVALFWALCVVCIVNMARYFSSLRALLVVLNRNCDPLLYQYVDGGGFFTSH GQPNKQV RLVWYIYAQRYRDHHDDEFIRRCERVRRQFILTSALCGLVVSLIALMIWH
Source	in vitro E.coli expression system
Target Names	uspB
Protein Names	Recommended name: Universal stress protein B
Expression Region	1-111
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.