



Recombinant Schizosaccharomyces pombe Putative uncharacterized protein C12C2.14c (SPBC12C2.14c)

Product Code	CSB-CF861204SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9C0X4
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	MEELQYNFKKRRKTHNGISRFRSALPLTIVYTIWSTFGSPCSGDQRVTLTSITSLRKVQ DRRESEKKVKGKGREEYRRYYFLLFYVSFPHIFLGLFFFIDKKILPFQSV
Source	in vitro E.coli expression system
Target Names	SPBC12C2.14c
Protein Names	Recommended name: Putative uncharacterized protein C12C2.14c
Expression Region	1-111
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.