



Recombinant Chicken Succinate dehydrogenase [ubiquinone] cytochrome b small subunit, mitochondrial (SDHD)

Product Code	CSB-CF719989CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5ZIS0
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	SHGGAPQGHGSSKAASLHWTSERAVSALLLGLLPAAYLYPGPAVDYSLAAALT LHGHWGL GQVITDYVHGDTPIKVANTGLYVLSAITFTGLCYFNYYDVGICKAVAMLWSI
Source	in vitro E.coli expression system
Target Names	SDHD
Protein Names	Recommended name: Succinate dehydrogenase [ubiquinone] cytochrome b small subunit, mitochondrial Short name= CybS Alternative name(s): Succinate dehydrogenase complex subunit D Succinate-ubiquinone oxidoreductase cytochrome b small subu
Expression Region	46-157
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.