



Recombinant Mesocricetus auratus Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit DAD1 (DAD1)

Product Code	CSB-CF006487MRG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61806
Product Type	Transmembrane Protein
Immunogen Species	Mesocricetus auratus (Golden hamster)
Sequence	SASVSVISRFLEEYLSSTPQRLKLLDAYLLYILLTGALQFGYCLLVGTFFNSFL SGFI SCVGSFILAVCLRIQINPQNKADFQGISPERAFADFLFASTILHLVVMNFVG
Source	in vitro E.coli expression system
Target Names	DAD1
Protein Names	Recommended name: Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit DAD1 Short name= Oligosaccharyl transferase subunit DAD1 EC= 2.4.1.119 Alternative name(s): Defender against cell death 1 Short name= DAD
Expression Region	2-113
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.