



Recombinant Enterobacteria phage f1 Head virion protein G6P (VI)

Product Code	CSB-CF304001ECH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P69531
Product Type	Transmembrane Protein
Immunogen Species	Enterobacteria phage f1 (Bacteriophage f1)
Sequence	MPVLLGIPLLLRFLGFLLVTLFGYLLTFLKKGFGKIAIAISLFLALIIGLNSILVGYLS D ISAQLPSDFVQGVQLILPSNALPCFYVILSVKAAIFFDVKQKIVSYLDWDK
Source	in vitro E.coli expression system
Target Names	VI
Protein Names	Recommended name: Head virion protein G6P Alternative name(s): Coat protein D G6P
Expression Region	1-112
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.