



Recombinant Na (+)/H (+) antiporter subunit C1

Product Code	CSB-CF350263FKZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P60682
Product Type	Transmembrane Protein
Immunogen Species	Staphylococcus aureus
Sequence	MEIIMIFVSGILTAISVYLVLKSLIRIVMGTTLLTHAANLFLITMGGLKHGTVPIYE AN VKSYVDPIPQALILTAIVIAFATTAFFLVLAFRITYKELGTDNVESMKGVPEDD
Source	in vitro E.coli expression system
Target Names	mnhC1
Protein Names	Recommended name: Na(+)/H(+) antiporter subunit C1 Alternative name(s): Mnh complex subunit C1 Mrp complex subunit C1
Expression Region	1-113
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.