



Recombinant Enterobacteria phage If1 Head virion protein G6P (VI)

Product Code	CSB-CF525685ELA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O80298
Product Type	Transmembrane Protein
Immunogen Species	Enterobacteria phage If1 (Bacteriophage If1)
Sequence	MPVFLGLPVLARFIGWLAGALIAYVAKFFTLGIARIALAISLFLGLIIGLNGLLVSYLSD LTSVLPPEIASAVSYVVPANAAPCLYAIFSLKAAVFIFDVKDRIIGYLDWNKS
Source	in vitro E.coli expression system
Target Names	VI
Protein Names	Recommended name: Head virion protein G6P Alternative name(s): Coat protein D G6P
Expression Region	1-113
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.