



Recombinant Araneus ventricosus Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit DAD1 (DAD1)

Product Code	CSB-CF818203AHAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8I7Z2
Product Type	Transmembrane Protein
Immunogen Species	Araneus ventricosus (Orbweaver spider) (Epeira ventricosa)
Sequence	MGSSAFEVLTFFLKDYKANTPQKLKIIDAYLLYILLTGINQFLYCCLVGTFFPFNSF LSGF ISCVASFVLGVCLRLQVNPQNSSNFCGIPPERAFADFIFAHVVLHLVVMNFIG
Source	in vitro E.coli expression system
Target Names	DAD1
Protein Names	Recommended name: Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit DAD1 Short name= Oligosaccharyl transferase subunit DAD1 EC= 2.4.1.119 Alternative name(s): Defender against apoptotic cell death 1 Defender
Expression Region	1-113
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.