



Recombinant *Betula pendula* Defender against cell death 1 (DAD1)

Product Code	CSB-CF006487BSS
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9M3T9
Product Type	Transmembrane Protein
Immunogen Species	<i>Betula pendula</i> (European white birch) (<i>Betula verrucosa</i>)
Sequence	MARSTSKDAQALIQSLRSAYAATPSNLKIIDLTVVFAVFTALIQVVYMAIVGSFPF NAFL SGLLSCIGTAVLAVCLRIQVNKENKEFKDLAPERAFADFVLCNLVHLVIMNFLG
Source	in vitro E.coli expression system
Target Names	DAD1
Protein Names	Recommended name: Defender against cell death 1 Short name= DAD-1
Expression Region	1-115
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.