



Recombinant Citrus unshiu Defender against cell death 1 (DAD1)

Product Code	CSB-CF006487DTR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9ZWQ7
Product Type	Transmembrane Protein
Immunogen Species	Citrus unshiu (Satsuma mandarin) (Citrus nobilis var. unshiu)
Sequence	MARSTGKDAQALFHSLRSAYAATPTTLKIIDLYVGFVFTALIQVVYMAIVGSFP FNSFL SGVLSCVGTAVLAVCLRIQVNKDNKEFKDLPPERAFADFVLCNLVLHLVIMNFL G
Source	in vitro E.coli expression system
Target Names	DAD1
Protein Names	Recommended name: Defender against cell death 1 Short name= DAD-1 Alternative name(s): CitDAD-1-1
Expression Region	1-115
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.