



Recombinant Succinate dehydrogenase hydrophobic membrane anchor subunit (sdhD)

Product Code	CSB-CF841277EOD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8X9A9
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli O157:H7
Sequence	MVSNASALGRNGVHDFILVRATAIVLTLTYIIYMGFFATSGELTYEVWIGFCASA FTKVF TLLALFSILIHAWIGMWQVLTDYVKPLALRLMLQLVIVVALVYVIYGFVVVWGV
Source	in vitro E.coli expression system
Target Names	sdhD
Protein Names	Recommended name: Succinate dehydrogenase hydrophobic membrane anchor subunit
Expression Region	1-115
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.