



Recombinant *Xenopus tropicalis* Transmembrane protein 218 (tmem218)

Product Code	CSB-CF023812XBF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B0JYX6
Product Type	Transmembrane Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Sequence	MATTILGVGPGVFIIAVIWMVTLMLTVLLCRASGKARFWTVVVFTLALITTLILVF FPRA SQTPAPAKEMQIVDTFFIGRYFLISIMSVIFLGCLFFVFVYHILEPVYAKPIGIH
Source	in vitro E.coli expression system
Target Names	tmem218
Protein Names	Recommended name: Transmembrane protein 218
Expression Region	1-115
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.