



Recombinant Enterobacteria phage IKe Head virion protein G6P (VI)

Product Code	CSB-CF361120ECQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P03674
Product Type	Transmembrane Protein
Immunogen Species	Enterobacteria phage IKe (Bacteriophage IKe)
Sequence	MPALLGIPALIRFIMGLVPIAIGYFAKFLGMIITRNGLMASALIGAILSVVVSFSIQLL GD ALSSSMGGISADFGNLMSSVLPDGTTCITVIITTRIAVFVFDIKDRLLGIANKVI
Source	in vitro E.coli expression system
Target Names	VI
Protein Names	Recommended name: Head virion protein G6P Alternative name(s): Coat protein D G6P
Expression Region	1-116
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.